

BRIDGE STAR
Safe Turning and

Repositioning System

**BridgeSTAR**<sup>™</sup> **System,** is a single patient use safe turning and repositioning system used by caregivers to safely manage a variety of frequent and strenuous patient handling tasks. The Star system will promote positive outcomes related to injury prevention and hospital-acquired pressure ulcers by creating a safe and efficient environment for both patients and staff.

The BridgeSTAR<sup>™</sup> is available in both air and sheet systems. The system can include foam wedges and highly breathable absorbent protection pads. Both applications are available in a variety of sizes and options. To make ordering and inventory management simple, all items of the system can be purchased separately.

## Why BridgeSTAR™ Systems?

- ✓ Injury Prevention Prevents injuries by significantly reducing the effort and forces required to perform turning and repositioning tasks
- ✓ Wound Care Incorporates breathable materials allowing the system to be left underneath the patient for high staff compliance when following turning and repositioning protocols. Eliminates sheering and friction
- ✓ **Sizing** Available in multiple sizes to meet the needs of your patients and the environment to achieve the highest compliance
- ✓ **Cost** It's a small price to pay to prevent injuries and health-acquired pressure ulcers. When you compare it to other costly solutions, it makes sense for your budget
- ✓ Versatile Effective in general Medsurge, Emergency, ICU's, and many other Departments and settings



T (800) 815-6615
E info@pelstarllc.com
W www.bridgehcusa.com



# Safe Turning and Repositioning System



#### BridgeSTAR™ Repositioning Sheets

- 800 lb weight capacity in standard sizes, 1000 lb weight capacity in x-large size
- Reduces force required to reposition patients
- Multiple handles to make turning simple

MODEL	DESCRIPTION	
PS-3035	BridgeSTAR™ Sheet Single Patient Use, 30" x 35"	12/cs
PS-3035VS	BridgeSTAR™ Sheet Single Patient Use, 30" x 35", w/ Velcro Straps	12/cs
PS-3555	BridgeSTAR™ Sheet Single Patient Use, 35" x 55"	12/cs
PS-3555VS	BridgeSTAR™ Sheet Single Patient Use, 35" x 55", w/ Velcro Straps	12/cs
PS-4278VS	BridgeSTAR™ Sheet Single Patient Use X-Large, 42"x 78", w/ Velcro Straps	12/cs

### BridgeSTAR™ Turning/ Repositioning Wedges

- Supportive wedges to reduce the risk of pressure ulcers
- Non-slip base ensures easy placement
- Directional arrows showing proper placement and position



MODEL	DESCRIPTION	
PS-W16115	BridgeSTAR™ Turning/ Repositioning Wedge, 25°, 16"x 11"x 5"	12/cs
PS-W16127	BridgeSTAR™ Turning/ Repositioning Wedge, 30°, 16"x 12"x 7"	12/cs
PS-W24127	BridgeSTAR™ Turning/ Repositioning Wedge, 30°, 24"x 12"x 7"	12/cs
PS-W34127	BridgeSTAR™ Turning/ Repositioning Wedge, 30°, 34"x 12"x 7"	6/cs
PS-W1788	BridgeSTAR™ Turning/ Repositioning Wedge, 45°, 17"x 8"x 8"	12/cs

#### BreatheDRY™ Absorbent Pads



- Soft, non-woven air-permeable top layer for patient comfort
- Manages climate control and preserves skin integrity



MODEL	DESCRIPTION	
BH-BD8040	BreatheDRY <sup>™</sup> Protect Full, 39" x 78", BridgeAIR <sup>™</sup> Absorbent Protection Pad	20/cs
BH-BD3135	BreatheDRY™ Standard 31" x 35", BridgeAIR™ Absorbent Protection Pad	60/cs
BH-BD4839	BreatheDRY™ Plus 48"x 39", BridgeAIR™ Absorbent Protection Pad	30/cs

# BridgeAIR<sup>™</sup> Mattresses



- Significant reduction in the forces required to laterally transfer and reposition patients
- Breathable fabric promotes a healthy microclimate for skin care



MODEL	DESCRIPTION	
PS34SPUL	BridgeAIR™ Mattress Single Patient Use Full, 34" x 78"	10/cs
PS39SPUL	BridgeAIR™ Mattress Single Patient Use Full, 39" x 78"	10/cs
PS50SPUL	BridgeAIR™ Mattress Single Patient Use Full, 50" x 78"	5/cs

