



BRIDGEAIR™

Single Patient Use
Lateral Transfer Solution

BridgeAIR™, an innovative air-assisted breathable lateral transfer and repositioning mattress, is designed to prioritize the well-being of caregivers, patients, and healthcare organizations. With its exceptional features, this mattress revolutionizes patient transfer by significantly reducing the physical effort required, offers advanced microclimate testing for optimal comfort, and provides value-driven purchasing programs that set it apart as the premier choice.

Design

- Quality-focused design and construction
- Rigorous clinical testing conducted with a leading independent laboratory
- Clinicians can trust in the reliability and effectiveness of BridgeAIR™ for patient care



Safety

- Reduce lateral transfer force by 80%-90%*
- Minimize staff and patient injuries related to lateral transfers and repositioning
- Easy to use; drives compliance



Microclimate

- Advanced materials allow for optimal breathability
- High-wicking properties maximize water vapor transmission rates
- Safe to leave underneath patient when not in use



Education and Training – B-SAFE Program – Contact us for more information.

*Data available upon request



PELSTAR LLC

T (800) 815-6615
E info@pelstarllc.com
W www.bridgehcusa.com

9500 W 55th St
McCook, IL 60525

PRODUCT SPECIFICATIONS

Width
34", 39", or 50"

Length
Full Size 78" or Half Size 48"

Material
Nylon, Latex free

Weight Limit
1200 lb/ 545 kg

Radiolucent and MRI Compatible

X-Ray: BridgeAIR™ Mattress Single Patient Use is radiolucent

MRI: BridgeAIR™ Mattress Single Patient Use is MRI compatible

(air supply must remain outside of MRI environment with extended supply hose)

BridgeAIR™ Air Supply

- Air-Assisted Transfer Device branded air supply pump manufactured to specification
- CE/UL 60601-1/CSA 601.1 M90 Certified by laboratory
- Hospital Grade Plug, 6'/25' Air Supply Hose

Now Available with Optional Bed Connecting Straps!



BridgeAIR™ Air Supply and BridgeAIR™ Transport Cart



BridgeAIR™ Mattress
Single Patient Use Half, 34" x 48"



BridgeAIR™ Mattress
Single Patient Use Full, 34" x 78"

ORDERING INFORMATION

MODEL	DESCRIPTION	
PS34SPUL	BridgeAIR™ Mattress Single Patient Use Full, 34" x 78"	10/cs
PS39SPUL	BridgeAIR™ Mattress Single Patient Use Full, 39" x 78"	10/cs
PS50SPUL	BridgeAIR™ Mattress Single Patient Use Full, 50" x 78"	5/cs
PS34SPUS	BridgeAIR™ Mattress Single Patient Use Half, 34" x 48"	10/cs
PS39SPUS	BridgeAIR™ Mattress Single Patient Use Half, 39" x 48"	10/cs
PS34SPULVS	BridgeAIR™ Connect Mattress Single Patient Use Full, 34" x 78" with Velcro Straps	10/cs
PS39SPULVS	BridgeAIR™ Connect Mattress Single Patient Use Full, 39" x 78" with Velcro Straps	10/cs
PS50SPULVS	BridgeAIR™ Connect Mattress Single Patient Use Full, 50" x 78" with Velcro Straps	5/cs
BHUT120	BridgeAIR™ Air Supply	1/ea
BHUARCRT	BridgeAIR™ Transport Cart	1/ea
BH-BD8040	BreatheDRY™ Protect Full, 39" x 78", BridgeAIR™ Absorbent Protection Pad	20/cs
BH-BD3135	BreatheDRY™ Standard 31" x 35", BridgeAIR™ Absorbent Protection Pad	60/cs
BH-BD3756	BreatheDRY™ Standard 37" x 56", BridgeAIR™ Absorbent Protection Pad	60/cs
BH-BD4839	BreatheDRY™ Plus 48" x 39", BridgeAIR™ Absorbent Protection Pad	30/cs

BridgeAIR™ Mattresses are also Available in Reusable Sizes



PELSTAR LLC

T (800) 815-6615
E info@pelstarllc.com
W www.bridgehcusa.com

9500 W 55th St
McCook, IL 60525